

Product information

MediaStrobo Version 3.5

New and improved functions



Main functions

- Microsoft Windows™ VISTA compatibility
- Compare-function for videos
- Compare-function for pictures
- Multiple-split display during recording and playing
- Recorder with display of current individual image recording next to live image
- Individual detail adjustment for image viewing
- Saving of patient-related document and video collections as data set
- Extended picture import dialogue alternative to recording

Ergonomics & subordinate functions

- Global or patient-related automatic generation of documents
- Direct control of record/stop during video recording without confirmation
- Automatic removal of black margins at image print-outs
- Fluent image display during rewind
- Simplified procedure of demand and input of the licence information
- GDT export with unambiguous data set assignation for direct selection
- Pre-selection of data sets to be compared in session mode/data set mode
- Manual adding of data sets to be compared from the entire player database
- Manual merging of split displays
- Automatic export of split displays in document collection
- Conversion of the online help system to HTML help system

Safety

- Operation without administrator authorisation (Windows™ VISTA)
- Automatic check for availability of a remote control unit
- Automatic configuration of audio and video devices
- Fault-tolerant intelligent input of date information

... and many other improved details such as:

- Optimised aspect ratio (width to height) in image display
- Complete deletion of the disk image before video export
- Optimised progressive image display in split view
- Automatic programme change upon GDT control