

Product Information

MediaStrobo Version 3.6

New and improved functions



Functionality & Ergonomics

- Screen keyboard for text input during the examination
- Automatic position adjustment of screen keyboard
- DICOM image export
- Double footswitch for separate recording of videos and individual images
- Managing of specialised databases for findings and diagnoses
- Extended support of USB cameras
- Individual selection of the image format when using video devices
- Uninterrupted video display during recording start/stop change
- Extended messages for user information
- Complete process monitoring upon archiving
- Complete export of databases in text file format for statistics
- Support of the „snap shot“ key of many USB video converters for recording control

Safety

- Microsoft Windows VISTA Service Pack 1 compatibility
- Automatic error detection and correction of database errors
- Improved network operation safety
- Database clearing function for the removal of deleted datasets

... and many other improved details such as:

- Audio recording also for USB video devices
- Minimal screen updating for a „smooth“ programme display
- Extended machine-analysable system records